S-PREMIUM

CABINA AUDIOMÉTRICA
AUDIOMETRIC TEST BOOTH

MANUAL DEL USUARIO
USER'S MANUAL

501-900-MUM · Rev. 1.02 · 2021-04

Manual Cabina Audiométrica **S-PREMIUM**

Revisión: 501-900-MUM Rev. 1.02

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SIBEL S.A.U.

Rosellón 500 bajos, 08026 BARCELONA (Spain)

Ventas Nacionales: Tel. 93 436 00 08 e-mail: comercial@sibelmed.com

International Sales: Tel. +34 93 436 00 07 e-mail: export@sibelmed.com

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Revision: 501-900-MUM Rev 1.00

All rights reserved

This manual can be purchased through the After Sales Service.

SIBEL S.A.U.

Rosellón 500 bajos, 08026 BARCELONA (Spain)

National Sales: Tel. +34 93 436 00 08 e-mail: comercial@sibelmed.com

International Sales: Tel. +34 93 436 00 07 e-mail: export@sibelmed.com

Technical Service: Tel. +34 93 433 54 50 e-mail: sat@sibelmed.com

Fax: +34 93 436 16 11, Web: www.sibelmed.com

ADVICE ABOUT THE COPYRIGHT

No part of this publication may be reproduced, transmitted, transcribed, stored in a back-up system or translated into any language or computer language in any form or by any means, electronic, mechanical, optical, chemical or manual without the express written consent from SIBEL S.A.U.

DISCLAIMER

SIBEL S.A.U. is responsible for the security, reliability and performance of this product only if:

- The place where the system is installed meets the requirements of electric installation IEC and other applicable regulations.
- Repairs, revisions or modifications, within and outside the warranty period are carried out by technical personnel of SIBEL S.A.U.
- The product is used by qualified staff in accordance with the recommendations of this User's Manual

**PRODUCT IN COMPLIANCE WITH MEDICAL
DEVICE REGULATION (EU) 2017/745.
(CLASS I)**

Revised
Date: 2021-04
Technical Director

Approved
Date: 2021-04
Commercial Director

INDEX

1. SAFETY	4
1.1. INTENDED USE	4
1.2. INDICATIONS FOR USE	4
1.3. WARNINGS.....	5
1.4. DISPOSAL OF ELECTRICAL AND ELECTRONIC EQUIPMENT BY DOMESTIC USERS IN THE EUROPEAN UNION	5
2. INSTALLATION AND INSTRUCTIONS FOR USE	6
2.1. MODELS AND OPTIONS	6
2.2. PAKING LIST.....	6
2.3. INSTALLATION AND SERVICING	6
2.4. INSTRUCTIONS FOR USE.....	7
2.5. CLEANING INSTRUCTIONS	7
3. TECHNICAL SPECIFICATIONS	8
4. APPLICABLE DIRECTIVES AND STANDARDS	10
5. SYMBOLS	11
6. PREVENTIVE MAINTENANCE.....	12

1. SAFETY

The **S-PREMIUM** booth has been designed and manufactured according to the Quality manual of SIBEL S.A.U., which agrees with the quality norms of EN 13485 e ISO 9001 and European Medical Device Directive 93/42/EEC. According to this directive, this is a Class I product.

The materials and finishing of the booth are made with the highest quality and has been tested by SounLab laboratories.

S-PREMIUM has been designed to guarantee the maximum safety. The instruction for use must be read before its use. Failure to do so may result in injuries to the user or patient and damage to the product and/or accessories. However, accidents may occur in the transport or storage, so it is recommended an initial check up of the booth and its accessories before the installation.

1.1. INTENDED USE

Test booth for the attenuation of external ambient noise intended for audiometric tests aiming at the detection, diagnosis, monitoring or treatment of hearing loss. It is intended to be used in combination with an audiometer. Not intended for outdoor use.

1.2. INDICATIONS FOR USE

The audiometric test booth **SIBELMED S-PREMIUM** is an individual chamber that attenuates the sound waves. This attenuation characteristics allows making audiometric tests according to standards ANSI S3.1 and ISO 8553-1. These standards specify that the audiometric tests must be performed in conditions where the external noise does not mask the results.

The booth can be used by patients older than 3 years old that are able to follow the instructions of the medical staff and have an appropriate mental state that allows being closed in a reduced enclosure during the test. The booth can be used by people in wheelchair with the optional ramp.

The medical staff must have the training in audiometric tests, know the functioning of the booth before its use and instruct the patient about the use of the booth.

1.3. WARNINGS

IF THE PACKAGING IS DAMAGED, CONTACT THE TRANSPORT AGENCY IMMEDIATELY AND YOUR DISTRIBUTOR BEFORE INSTALLING THE BOOTH. DO NOT DISPOSE THE PACKAGING UNTIL THE CORRECT FUNCTIONING OF THE PRODUCT HAS BEEN THOROUGHLY CHECKED.

Ensure the door is sealed before each audiometric test.

Use the fairlead, placed next to the table, to connect the audiometer and the accessories. Once the wires have been leaded by the fairlead, it must be screwed to avoid attenuation loss.

The lateral gutters can be used to connect the speakers and make the free field tests. The final connection will be done depending on each audiometer.

Report any serious incident to the manufacturer and competent authority.

1.4. DISPOSAL OF ELECTRICAL AND ELECTRONIC EQUIPMENT BY DOMESTIC USERS IN THE EUROPEAN UNION

The **S-PREMIUM** booth must be disposed properly; it may need to be recycled in accordance with the statutory requirements in your country.

- Materials according RoHS Directive: the product and its accessories are RoHS compliant.
- Electrical and electronic waste according to WEEE Directive (2002/96/CE): the user is the responsible of disassembling the electrical and/or electronic components of the booth and delivering them in a recollecting point of electric and electronic devices.

Information on proper disposal is available from SIBEL's After Sales Service or your dealer.

2. INSTALLATION AND INSTRUCTIONS FOR USE

2.1. MODELS AND OPTIONS

The S-PREMIUM Booth has the following models:

Reference	Model	Dimensions
09041	A	129 x 129 x 212 cm
09042	B	138 x 106 x 222 cm
09043	C	198 x 168 x 222 cm
09044	D	198 x 198 x 222 cm
09045	E	228 x 198 x 222 cm
09046	F	228 x 228 x 222 cm
09047	G	258 x 228 x 222 cm
09048	H	258 x 258 x 222 cm
09049	I	288 x 258 x 222 cm
09050	J	288 x 288 x 222 cm
09051	K	Special dimensions

2.2. PAKING LIST


The audiometric booth is composed by the next elements and accessories:

CODE	QUANT.	DESCRIPTION
09041-09051	1	BOOTH S- PREMIUM MODEL A-K
09165	1	USER'S MANUAL
According to the model	1	ASSEMBLY GUIDE
09052	1	FORCED VENTILATION (Optional)
09053	1	ACCES RAMP (Optional)
09054	1	WASHABLE FLOOR (Optional)
09109	1	DOOR OPENING TOWARDS LEFT SIDE (Optional)
----	1	WINDOW AT RIGHT SIDE OF DOOR (Optional)

2.3. INSTALLATION AND SERVICING

Consult the Assembly Quick Guide to assemble and put into service the booth.

2.4. INSTRUCTIONS FOR USE

The booth has a door with lock (1), a window that allows the mutual vision between patient and user (2), illumination activated by a movement sensor (3), an external table for placing the audiometer (4) and a fairlead for the connections of the accessories of the audiometer (5). Also, the booth includes absorbent panels to avoid the coincidence effect (6) and holders below the inner table (7) to hold the headphone, patient's pushbutton etc, when not in use.

2.5. CLEANING INSTRUCTIONS

The external surface and the inner walls can be cleaned with soap and water, drying the residual moisture. The carpeted floor can be cleaned by a conventional vacuum cleaner, the plastic parts and the window can be cleaned with water and soap. The absorbent panels can be cleaned with neutral soap and water. Do not use abrasives substances or solvents that can damage the finishing of the booth.

3. TECHNICAL SPECIFICATIONS

	MODEL A	MODEL B	MODEL C	MODEL D
Inner dimensions (cm)	113 x 113 x 195	122 x 90 x 205	164 x 134 x 210	164 x 164 x 210
External dimensions (cm)	129 x 129 x 212	138 x 106 x 222	198 x 168 x 222	198 x 198 x 222
Window dimensions (cm)	82 x 60	60 x 60	41 x 60	60 x 60
Inner table (cm)	113 x 11,5 x 3	90 x 11,5x3	70 x 11,5 x 3	70 x 11,5 x 3
Door entrance width	71cm	65cm	80cm	80cm


	MODEL E	MODEL F	MODEL G	MODEL H
Inner dimensions (cm)	194 x 164 x 210	194 x 194 x 210	224 x 194 x 210	224 x 224 x 210
External dimensions (cm)	228 x 198 x 222	228 x 228 x 222	258 x 228 x 222	258 x 258 x 222
Window dimensions (cm)	60 x 60	60 x 60	60 x 60	60 x 60
Inner table (cm)	70 x 11,5 x 3	70 x 11,5 x 3	70 x 11,5 x 3	102 x 11,5 x 3
Door entrance width	80cm	80cm	80cm	80cm

	MODEL I	MODEL J
Inner dimensions (cm)	254 x 224 x 210	254 x 254 x 210
External dimensions (cm)	288 x 258 x 222	288 x 288 x 222
Window dimensions (cm)	60 x 60	60 x 60
Inner table (cm)	102 x 11,5 x 3	102 x 11,5 x 3
Door entrance width	80cm	80cm

Medical device classification	Classe I
Ambient use conditions	Temperature: 5 up to 40 °C Humidity: 5% - 85 % (Without condensation)
Transport and storage conditions	Temperature: -5 up to 45 °C Humidity: 5% - 65% (Without condensation)
Lifetime	7 years
Minimal Global Attenuation for white noise between 250 and 8000Hz	Depending on the model
Attenuation Global maximal for white noise	Depending on the model
Door	Glass
Lock points	4 hinges
Window type	Double glass with PVC frame
Base	Silent block
Outer table	Fixed and removable
Table dimensions	113 x 45 x 3 cm.
External and internal finishing	Varnished
Absorbent panels	Yes
Floor	Carpet
Ventilation	Free by laberinth in the roof
Illumination (*)	LED light - 220V
Free field pre-installation	Yes
Forced ventilation	(Optional)
Acces Ramp	74 x 74 x 10 cm (Optional except Model B)
Floor	Washable (Optional)

(*) 110V for other countries.

Frequency <i>f</i> Hz	D'p (octave third) dB
100	18,5
125	24,2
160	16,6
200	22,5
250	31,4
315	23,5
400	23,6
500	31,9
630	33,6
800	33,7
1000	34,4
1250	34,9
1600	34,6
2000	33,4
2500	35,7
3150	36,9
4000	39,3
5000	41,2
6300	42,1
8000	43,6
10000	43,4

Test performed according to EN ISO 11957:2009

4. APPLICABLE DIRECTIVES AND STANDARDS

- Medical device regulation (EU) 2017/745
- RoHs Directive 2011/65/UE
- Directive 2012/19/UE on waste of electrical and electronic equipment (WEEE) (RD 110/2015)
- Audiometry: Allows the fulfillment of ISO 8253-1 and ANSI S3.1 depending on surrounding conditions
- Quality: Sibel's Quality System (EN ISO 9001:2015 and EN ISO 13485:2016+AC:2018)
- Safety: EN 60601-1:2006+AC:2010 + A1:2013, EN 60601-1-6:2010 + A1:2015; EN 62366:2008 + A1:2015
- Risk management: EN ISO 14971:2019
- Information provided by the manufacturer: EN 1041:2008 + A1:2013
- Symbols used: EN 15223-1:2016
- Acoustic insulation: EN ISO 11957:2009

5. SYMBOLS

MANUFACTURER

DATE OF MANUFACTURE

CATALOGUE NUMBER

SERIAL NUMBER

MEDICAL DEVICE

CE MARKING

CAUTION, ADDITIONAL WARNINGS IN ACOMPANYING DOCUMENTS

CONSULT INSTRUCTIONS FOR US

WARNING, RISK IDENTIFIED

WASTE OF ELECTRICAL AND ELECTRONIC EQUIPMENT

TEMPERATURE LIMITS

HUMIDITY LIMITATION

KEEP DRY

FRAGILE, HANDLE WITH CARE

THIS WAY UP, PACKAGING

6. PREVENTIVE MAINTENANCE

The user should check the aspect of the booth periodically, the internal illumination, the locking door, confirming that there is no part or element damaged.

If there is any abnormality that does not allow the normal use of the booth, contact with the After Sales Service of Sibel S.A.U., specifying the problem as accurately as possible.