# **S3000**

**JERINGA DE CALIBRACIÓN**

**CALIBRATION SYRINGE**

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# Index


<b>1. Indications for Use .....</b>	<b>14</b>
<b>Syringe diagram .....</b>	<b>15</b>
<b>Intended Use and Indications for Use .....</b>	<b>15</b>
<b>Directions For Use .....</b>	<b>16</b>
<b>Reprocessing Instructions .....</b>	<b>17</b>
<b>Scope.....</b>	<b>17</b>
<b>Product classification.....</b>	<b>17</b>
<b>Decontamination .....</b>	<b>17</b>
<b>Cleaning Agents/Supplies Calibration Syringes ....</b>	<b>18</b>
<b>Cleaning/Maintenance Methods .....</b>	<b>18</b>
<b>Step 1. Inspection.....</b>	<b>18</b>
<b>Step 2. Manual cleaning (waste disposal) .....</b>	<b>18</b>
<b>Step 3. Functional test.....</b>	<b>18</b>
<b>Step 4. Recalibration.....</b>	<b>18</b>
<b>Step 5. Storage .....</b>	<b>18</b>
<b>Symbols .....</b>	<b>19</b>
<b>Useful life.....</b>	<b>19</b>
<b>Warning .....</b>	<b>20</b>
<b>Safety Information .....</b>	<b>21</b>


# 1. INDICATIONS FOR USE

## Syringe diagram

### Option 1

### Option 2

## Intended Use and Indications for Use

SIBELMED Volume Calibration Syringes are designed for testing accuracy and if it's necessary, to calibrate spirometers. These Calibration Syringes are manually operated by actuating the piston rod in and out of the cylinder assembly. The outlet port of the Syringe assembly is attached to the spirometer transducer. The Calibration Syringe, when actuated, displaces known volumes of air into the flow or volume measurement device. The operator compares the calibration syringe volume with the results obtained from integrating the flow or reading the displaced volume of the measurement device to determine its accuracy. Calibration Syringes are used in pulmonary function test laboratories (both hospitals and clinics), bio-medical engineering calibration laboratories, private physician offices, universities,

research centers, and by manufacturers of respiratory systems. SIBELMED syringes are used by trained staff to perform the calibration procedures on the spirometers or other equipment under test.

Warning: It's very important that the operator performing the calibration with the syringe is trained. Otherwise measurement errors can be introduced in the spirometer.

## Directions For Use

For detailed information contact SIBEL S.A.U.

The SIBELMED syringe is typically used to calibrate spirometers and has a volumen of 3 liters. Its accuracy is  $\pm 0.5\%$ .

Standards and recomendations:

- Medical device regulation (EU) 2017/745: CE Mark Class I
- ATS / ERS 2005: Miller et al. Task force: Standardisation of Lung Function Testing. Standardisation of Spirometry. Eur Respir J 2005; 26:319-338
- EN 15223-1:2016: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied
  - Part 1: General requirements.
- EN 1041:2008 + A1:2013 Information supplied by the manufactures of Medical Devices.
- SIBEL system quality conforms with EN ISO 13485:2016+AC:2018, EN ISO 9001:2015 and EN ISO 14971:2019

Three liters are confirmed when the piston rod pulled all the way out and it stops at the zero position.

SIBELMED syringe has been calibrated sealing the locking collar. Therefore, it must never be moved under normal use.

Confirm that the connection between the syringe output port and the device under test is leak free.

Volume calibration of the spirometer is performed with the syringe by actuating the piston rod from the zero position to the output volume. The output volume is confirmed when the locking collar stops at zero position.

With the gray syringe model (Option 2), and for your calibration of Sibelmed spirometers, use the included calibration adapter (code 09497), as shown in the previous figure.

When the syringe is not in use, place the plug provided in the outlet port to keep the syringe internal volume clean and free of contaminants.

## Reprocessing Instructions

### Scope

This guidance is directed to personnel responsible for cleaning and maintaining SIBELMED calibration syringes.

### Product classification

The syringe is considered Class I according to Medical device regulation (EU) 2017/745.

SIBELMED calibration syringes are supplied clean, non-sterile. Calibration syringes fall into the non-critical device category based on potential risk of infection. Non-critical equipment should be kept visibly clean with soap and water or an intermediate level disinfectant. This device will not come into contact with patients therefore disinfection or sterilization is not required.

## Decontamination

SIBELMED calibration syringes are not required to be disassembled in the field for cleaning or disinfection. Without proper training and equipment the precise components of these syringes could be damaged in the disassembly/reassembly process. These calibration syringes should be returned to SIBEL S.A.U. once a year for maintenance and recalibration.

It is the responsibility of the user (health care personnel) for ensuring that: the cleaning methods recommended by SIBEL S.A.U. can be duplicated in their environment, that appropriate tools, and replacement parts are available and that instructions provided in this user's manual are followed correctly.

## Cleaning Agents/Supplies for Calibration Syringes

Mild detergent (neutral pH7) or soap with tap water or standard hospital disinfectant/cleaners are recommended for periodic cleaning (wiping down) of SIBELMED calibration syringes. Use clean soft cloths for wiping and drying.

Caution: Do not use cleaning solutions or cloths inside the cylinder area of the syringe.

## Cleaning/Maintenance Methods

### Step 1. Inspection

Inspect the syringe for damage at all stages of handling. If damage is detected on any components of the syringe it should be identified and documented. Complete the cleaning process and contact SIBEL S.A.U. for guidance on replacement components and after sales service.

### Step 2. Manual cleaning (waste disposal)

Prepare the detergent or soap solution. Typical concentration of detergent is 30g (one ounce) to 3,8 liters of water. With a clean soft cloth damp from the detergent solution wipe down the outside surfaces of the syringe. Follow with a dry soft cloth to remove any residue. When using hospital disinfectant cleaners follow the manufacturer's instructions. Dry thoroughly using a soft cloth or disposable paper towel.

### Step 3. Functional test

Leak test the calibration syringe by occluding the outlet port with one hand while pushing on the piston rod with the other. With your ear close to the piston rod end plate listen for any hissing sound indicating a leak.

### Step 4. Recalibration

Recalibration at SIBEL S.A.U. is recommended once per year. A calibration due date is labeled on each syringe.

### Step 5. Storage

Calibration syringes should be stored in a way that prevents contamination and damage between uses.

Storage conditions:

Temperature: -40 – 70 °C

Humidity: 0 – 95 %

1. Place the plastic plug into the syringe outlet port to prevent contamination of the internal volume of the syringe. Each syringe is shipped with a plastic plug in the outlet port.

2. Place the syringe assembly in a clean plastic bag and close the end.

When maintained properly and recalibrated on a regular basis (annual), SIBELMED calibration syringes should optimally function within specifications. Visual inspection of all the components of the assembly along with the functional tests as described will identify a defective syringe.

## Symbols

MANUFACTURER

(manufacture date, name and address of the manufacturer)

CE MARKING

MEDICAL DEVICE

CATALOGUE NUMBER

SERIAL NUMBER

TEMPERATURE LIMITS

## Useful life

The useful life of the syringe is 5 years. When maintained properly and recalibrated on a regular basis (annual), SIBELMED calibration syringes should optimally function within specifications. Visual inspection of all the components of the assembly along with the functional tests as described will identify a defective syringe.

## Warning

1. SIBELMED calibration syringes are not intended for direct or indirect connection to patients (or animals) for purposes of conducting diagnostic tests or for ventilation or generation of pressure-volume curves by insufflations of a patient's (or animal) lungs. There are no safety mechanisms built into the syringes and the inside of the syringes cannot be disinfected without damage.
2. Do not attempt to sterilize these calibration syringes.
3. Do not disassemble this product and avoid liquid or dirt to get inside of it. Contact SIBEL S.A.U., for information on spare parts, recalibration or after sales service.
4. Do not operate at temperatures outside 5 to 40 °C range.
5. Potential hazard of incorrect calibration of diagnostic flow, volume transducer or device being calibrated.
6. Potential body injury if the syringe rolls, falls from the table or is used in a wrong way (avoid putting the fingers in the pistol travel range when moving the piston)
7. Potential contamination hazard of syringe from hospital or other contaminants. Devices to be calibrated with a syringe must be disinfected in order to prevent the syringe being contaminated.
8. The disposal of this product must be carried out in conformity with CE 1272 /2008
9. Storage conditions must be observed.
10. If you find any damage to the packaging, please contact the shipping agent and your dealer immediately prior to installing or using the syringe. Do not get rid of the packaging until the proper operation of the syringe is verified.

Report any serious incident to the manufacturer and competent authority.

## **Safety Information**

To request technical or safety information on SIBEL S.A.U. calibration syringes contact SIBEL S.A.U. by calling +34 93 433 00 08 or by fax to +34 93 436 16 11. You can also contact directly the dealer in the area.

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